



# COVID-19 Cleaning Products



**HOW**  
Diluted to spray bottle with mixing unit at station

**WHERE**  
Food prep surfaces, mask/regulators and apparatus

**WHEN**  
Twice daily and after exposure



SARS Associated Coronavirus    547 PPM    10 Minutes    ¼ oz/gl



**HOW**  
Diluted to spray bottle with mixing unit at station

**WHERE**  
Work surfaces, non-porous items and apparatus

**WHEN**  
Twice daily and after exposure



SARS Associated Coronavirus    850 PPM    10 Minutes    2 oz/gl



**HOW**  
Use per product labeling instructions

**WHERE**  
All EMS equipment, other on-scene items and apparatus

**WHEN**  
Daily and after every call



**FOR USE AS A VIRUCIDE\* ON NON-INSTRUMENT SURFACES, INCLUDING HIV-1, HBV AND HCV APPLICATIONS:**  
Spray **CaviCide** directly onto pre-cleaned surface, thoroughly wetting area to be treated. Allow surface to remain visibly wet for 2 minutes at room temperature (68 °F/20°C).



**HOW**  
Make a 10% bleach solution

**WHERE**  
Floors, bathrooms and if Cavicide or 764 runs out

**WHEN**  
Twice daily and after exposure



Human Coronavirus    2500 PPM    5 Minutes    6 oz/gl